

**HDH**  
*Clinic Patients*  
**Specimen Collection Locations**

- Brock 1 clinic
- Jeanne Mance 5 clinic

**Requisition completed by:** \_\_\_\_\_

**Sample collected by:** \_\_\_\_\_

**Physician Name (REQUIRED):** \_\_\_\_\_

Location: \_\_\_\_\_

Date/Time of Collection: \_\_\_\_/\_\_\_\_/\_\_\_\_ / \_\_\_\_:\_\_\_\_

yyyy mm dd hhmm

Material Submitted:  Blood  Urine  Fluid: \_\_\_\_\_ (specify)

**Light Green (Li Heparin PST)**

- Albumin
- Alcohol (Ethanol)
- Alkaline Phosphatase (ALP)
- ALT
- BHCG, Total (HCG)
- Bilirubin, Total
- Calcium, Total
- Cholesterol, Total
- CK, Total
- Cortisol  AM  PM
- Creatinine (includes eGFR)
- CRP (High Sensitivity)
- Electrolytes (Na, K, Cl, TCO<sub>2</sub>)
- Estradiol
- Ferritin
- FSH
- Glucose, Fasting
- Glucose, Random

**Glucose Tolerance Tests\***

- 50 g Screen for Gestational Diabetes
  - 1 hour
- 75 g Confirmation of Gestational Diabetes
  - fasting  1 hour  2 hour
- 75 g test for Non-Pregnant Adults
  - fasting  2 hour

**\*Multiple collections require individual requisitions**

**Glucose Meter Split Sample**

- GM result \_\_\_\_\_ Operator ID: \_\_\_\_\_
- Haptoglobin
- Iron\*\*
- Iron % Saturation\*\* (Iron & Transferrin)
- LH (luteinizing hormone)
- Lipase
- Lipid Profile, Fasting (Chol/Trig/HDL/LDL)
- Magnesium
- Phosphate
- Prolactin
- Protein, Total
- Triglycerides
- Thyroxine, Free (FT4)
- TSH
- Transferrin
- Urea
- Uric Acid

\*\* A.M. Collection Recommended

**Gold (Gold Top SST)**

- AFP
- Beta 2 Microglobulin
- Beta Hydroxybutyrate (BOHB/Ketones)
- Calcium, Ionized (needs separate tube)
- CA 15-3 (breast)
- CA 19-9 (pancreas)
- CA 125 (ovarian)
- CEA
- C3
- C4
- DHEAS
- Folate, Serum (*protect from light*)
- IgE (allergy)
- Immunoglobulin, IgG
- Lactate Dehydrogenase (LDH)
- Osmolality
- Protein Electrophoresis (SPE) (needs separate tube)
- Progesterone
- PSA, Total
- PSA, Free
- Rheumatoid Factor
- SHBG
- Testosterone\*\*
- Testosterone\*\*, Bioavailable (calculated)
- Thyroid Peroxidase Antibodies
- Thyroglobulin
- Thyroglobulin Antibodies
- Vitamin B12

**Hep Syringe or Pediatric Capillary Tube on Ice Water**

- Blood Gases
  - FiO<sub>2</sub> \_\_\_\_\_
  - Arterial
  - Capillary
  - Venous

**Hep Syringe, Dark Green (Li Heparin, No Gel) or Pediatric Capillary Tube on Ice Water**

- Carboxyhemoglobin
- Methemoglobin
  
- Lactate

**Order of Draw**

- Blue (Na Citrate) \_\_\_\_\_
- Red (no gel) \_\_\_\_\_
- Gold (gel) \_\_\_\_\_
- Lt Green (Li Hep) \_\_\_\_\_
- Dark Green (Li Hep) \_\_\_\_\_
- Dark Green (Na Hep) \_\_\_\_\_
- Lavender (EDTA) \_\_\_\_\_
- Pink (EDTA) \_\_\_\_\_
- Royal Blue (K2 EDTA) \_\_\_\_\_
- Grey (Na Fluoride) \_\_\_\_\_
- Venous Blood Gas (Hep syringe) \_\_\_\_\_

**Lavender (EDTA)**

*Separate tube required for each test*

- ACTH – Requires 2 pre-chilled tubes (**on ice water**)
  - ADH – Requires 2 pre-chilled tubes (**on ice water**)
  - Ammonia (**on ice water**)
  - Cyclosporin \*
  - Folate, RBC (*protect from light*)
  - HbA1C
  - PTH (Parathyroid Hormone – Intact)
  - Sirolimus \*
  - Tacrolimus (FK506) \*
  - Troponin (includes CK – requires Li Heparin PST)
- \* Indicate date and time of last dose**
- \_\_\_\_/\_\_\_\_/\_\_\_\_ / \_\_\_\_:\_\_\_\_
- yyyy mm dd hhmm

**ROUTINE HEMATOLOGY (Lavender)**

- CBC (*Specials on back*)
- Differential
- ESR
- Retic Count
- Body Fluid Cell Count & Diff \_\_\_\_\_ (*collection site*)

**ROUTINE HEMOSTASIS (Specials on back)**

**(Blue 2.7 or 4.5 mL draw)**

- PT (includes INR)
- PTT
- PT 50/50 mix
- PTT 50/50 mix
- D-Dimer
- Fibrinogen
- Thrombin Time

**MEDICATION**

- Oral Anticoagulant:  No  Yes
- If yes, type: \_\_\_\_\_
  
- Heparin:  No  Yes
- If yes:  UFH or  LMWH
- Other: \_\_\_\_\_

